

Eclips+

STANDARD & LIGHT MANUAL WHEELCHAIRS

Lightweight aluminum wheelchair with a fixed back and reinforced cross. Equipped with Quickrelease rear wheels, foldable and removable armrests.



C20
blue

C25
diamond
black

C29
grey

FEATURES

A solid classic throughout the years



Optional (HEM2) hemiplegic kit



Optional (T30) push version




Ergonomic handles

Foldable and removable armrests

Rear wheels with Quickrelease system

Foldable footplates

Wide range of accessories and options.



590	300	390	220	430-530	400-460	390-460	420-440-460	1040-1080	14,5	130
620	300	420	220	430-530	400-460	390-460	420-440-460	1040-1080	14,5	130
640	300	440	220	430-530	400-460	390-460	420-440-460	1040-1080	15	130
660	300	460	220	430-530	400-460	390-460	420-440-460	1040-1080	15,8	130
680	300	480	220	430-530	400-460	390-460	420-440-460	1040-1080	15,8	130
700	300	500	220	430-530	400-460	390-460	420-440-460	1040-1080	15,9	130
750	300	550	220	430-530	400-460	390-460	420-440-460	1040-1080	16	130



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.