

Eco-Light II

Rollator

The **Eco-Light II** is a foldable 4-wheel rollator. It is equipped with a removable basket, padded seat and crutch holder. A foldable bar provides extra support when seated. Complete with parking brakes and height adjustable handles. Foldable and ready for transport in seconds.

total width	615 mm
seat depth	305 mm
seat width	355 mm
seat height	560 mm
total height	800 - 930 mm
total length	695 mm
weight	7,5 kg
max. users weight	120 kg
wheels (front/rear)	8 x 1 ^{1/4} (ø 200 mm)
colour(s)	grey

Measurement tolerances: ±15mm / 1,5kg.
We reserve the right to introduce technical changes.



N.V. VERMEIREN N.V. - Vermeirenplein 1-15 - B-2920 Kalmthout - Belgium
Tel.: +32 (0)3 620 20 20 - Fax: +32 (0)3 666 48 94 - info@vermeiren.be - website: www.vermeiren.com

Images for illustration purposes only. Certain features may be optional.

Opening up new horizons in mobility