

# 9062 - 9063

## Commode chair

**This classic commode chair with artificial leather upholstery delivers optimal comfort and hygiene.** Sturdy fixed armrests and a strong backrest ensure the premium quality of this commode chair. The backrest is enveloping for a better seating comfort, the padded upholstery feels soft to the touch and the distance between the seat and the bucket is minimal to ensure a better hygiene.



*Model 9063*  
Height adjustable legs



*Model 9062*





*Enveloping backrest for a better seating comfort*



*9063 - Height adjustable legs*



*Soft padding for a softer body contact*



*Minimal distance between the seat and bucket*

## Frame colour(s) Upholstery colour(s)



C86



## Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the **9062 and 9063** comply with the highest norms, laid down by certified test institutes in most European countries.

## Technical data

■ Total width	620 mm
■ Total height	9062 : 835 mm - 9063 : 755 - 885 mm
■ Total length	650 mm
■ Seat width	430 mm
■ Height armrests	215 mm
■ Seat height	9062 : 475 mm - 9063 : 470 - 560 mm

■ Height backrest	400 mm
■ Seat depth	450 mm
■ Weight	9062 : 8,35 kg - 9063 : 9 kg
■ Maximum user weight	120 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN Group N.V.

Vermeireplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer

**VERMEIREN**  
we care for you