

Forest 3 SU

Stand-up

The Forest 3 SU is a power wheelchair with stand-up possibility. Just like its brother, the Forest 3 but with additional functions.

Complete with a safety belt at chest level and a padded knee guard for complete comfort and safety.

The control unit allows the user to shift from the standing position to the sitting position to the 180° reclined position.

The small wheels on the front footplate provide additional stability when driving your Forest 3 SU.





Colours




Quality guaranteed

Regarding quality, the **Forest 3 SU** complies with the highest norms, laid down by certified test institutes in most European countries.

Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.



Technical data

Total width	602 / 670 mm
Total length	1170 mm
Total height	1100 - 1575 mm
Seat height	530 mm
Seat width	400 / 450 mm
Seat depth	420 - 470 mm

Batteries / Motors	2 x 12V / 80Ah - 220W
Speed	6 km/h
Driving range	50 km
weight	165 kg
Maxiumum user weight	130 kg
Additional information	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you