

# D200 30°

## Standard wheelchair with reclining backrest

D200 30° is a lightweight aluminum wheelchair with a 30° reclining backrest (adjustable without tools). A standard manual wheelchair perfect for a rental parc or personal use. D200 30° is a fully adaptable chair to suit any user or changes that are needed.



**VERMEIREN**  
we care for you

## Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.


All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.



Rear wheels with quick-release system



Removable swing-away leg rests (in and outwards)



Height adjustable (3 different heights)



Armrests adjustable in height and depth to suit specific needs


## Technical data

Total width	580/600/620/640/660/680/700/720 mm
Total length	1010 mm
Total height	880 - 940 mm
Folded width	300 mm
Seat width	380/400/420/440/460/480/500/520 mm
Height armrests	220 - 240 mm

## Quality guaranteed

Regarding quality, the **D200 30°** complies with the highest norms, laid down by certified test institutes in most European countries.

## Colours



Option B78 - Anti-tipping wheels

Option BZ8 - articulated legrest self-correcting in length



Seat height	470 - 530 mm
Height backrest	420 mm
Seat depth	430 mm
Weight	17,9 kg
Maximum user weight	130 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN Group

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer

**VERMEIREN**  
we care for you