

Stabilo

C U S H I O N S



Stabilo cushions

Introduction


The seat and backrest are made with vacuum technology. For the production, high quality materials such as flexible neopren and modified polyester granules are used so that the cushions will not deform, even with long term pressure.

First we fill the cushion with air. It is flexible and matches the anatomical shape of the body of the seated person, ensuring distribution of the pressure. Second, we remove the air from the interior cushions to permanently fix the shape.

Complete formation and shape of the surface of the cushion gives an excellent stability. The system creates a big comfort for people with a deformed spine. It is indispensable for people in need of anatomical support. Recommended for people with with spastic and anti-bedsore prophylactics.

Your guarantee of quality.

Regarding quality, our products comply with the highest norms, laid down by certified test institutes in most European countries.



N.V. Vermeiren N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

Tel.: +32 (0) 3 620 20 20
Fax: +32 (0) 3 666 48 94

www.vermeiren.com - info@vermeiren.be

For more information, contact your local dealer



VERMEIREN
Opening up new horizons in mobility

Subject to technical alterations

Information for the visually impaired is available upon request. - R.E.: N.V. Vermeiren N.V., Vermeirenplein 1/15 - 2920 Kalmthout - Belgium - 06/2014