





Thanks to years of experience in care wheelchairs, we accommodate to the youngest users and the caregivers.

The design of Inovys Junior is expertly thought out to suit individual needs. It brings comfort, versatility, with endless possibilities, on top of Junior friendly design.

Efficient postural management

For short or long term usage, postural management throughout the day is a must for our users. Either for pressure relief, blood circulation, or just for comfort. You can easily adjust the tilting system to suit your needs.

Gas spring tilting & inclination possibilities for more precision and more comfort for the Junior and caregiver.


Endless adjustment possibilities

Adjustment and configuration do not need to be set exclusively when comes the ordering process.

During the day, the use, or with user growth, adjustment possibilities are a must to fit updated needs and to keep providing comfort.

Optimal comfort solution

Inovys Junior version offers comfort pre-shaped seat and backrest cushions, available in many versions.



Individual needs configuration

Not only does the chair allows to vary position & configuration. You can also configure and order a Inovys Junior that fits the users' needs

best, thanks to the wide range of options and accessories.



Consult our website or order form to review all options.



Technical data

Total width	580 mm
Total height	990 - 1060 mm
Total length	1050 mm
Seat angle	0° - 21° mm
Seat width	350 mm
Height armrests	220 - 320 mm

Seat height	390 - 450 mm
Height backrest	550 - 750 mm
Seat depth	350 - 470 mm
Weight	36 kg
Maximum user weight	60 kg
Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com

For more information, contact your local dealer

