



* - used only with TURIOS LE + SE19

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

<h2>Armrest SE53 Turios</h2>				-
				Plan No. EP00488/J
Scale: 1:1	Drafted	Verified	Confirmed	
Date	12.02.2020	7.11.2025	7.11.2025	
Initials	DR	FM	MK	

Electrical components presented in 2992521/2992522/2992523/2992524/2992525

Material			
No.	Plan	Article	Quantity
1	-	2606454	1
2	Y003626	1476554	1
3	Y009175	1641654	3
4	Y009176	1641754	1
5	YP07839	2783154	2
6	YP08051	2605154	1
7A	YP08058	2605454 (R)	1
	YP08059	2605554 (L)	1
7B*	P00112681	2911954 (R)	1
	P00112683	2912054 (L)	1
8	-	-	-
9	YP09265	2642654	1
10	-	-	-
11	-	-	-
12	X005565	1912287	1
13	Armrest right	1913465	1
	Armrest left	1313464	1
14	XP00227	2990147	2
15	-	-	-
16	Washer M6	1000973	8
17*	M5x12	1001489	2
18	Tube Cap 20x20	1001612	2
19	Nut M10	1001672	1
20	Nut M5	1001679	4
21	Nut M6	1001680	6
22	M6x40	1003111	2
23	Star Knob	1003658	1
24	M5x20	1003907	4
25	M6x40	1006111	4
26	-	-	-
27	-	-	-
28	-	-	-
29	Nut M6	1008492	2
30	-	-	-

Material			
No.	Plan	Article	Quantity
31*	M4X12	2994126	2
32	-	-	-
	-	-	-
	-	-	-
	-	-	-
	-	-	-
	-	-	-
33	M5x16	1907493	2
34*	-	-	-
35	Spring plunger	1991432	1
36	-	-	-
37*	Washer M4 DIN127B	1998515	2
38	M6x45	2994058	2
39	Tube Cap 16x16	2994195	2
40	XP03080	2996656	1

Bill of Material 1913465/1913464			
No.	Plan	Article	Quantity
13a	X005567	1912286 (R)	1
	X005566	1912285 (L)	1
13b	X005517	2014999	1
13c	X005519	1911778	1
13d	X005561	2021699	1
13e	M5x20	1915452	2
13f	M6x16	2994067	1
13g	Nut M6	2994059	1
13h	Rivet Ø4	2995285	2
13i	M5x16	3992507	2
13j	Washer M5	2995288	2
13k	Sticker	T001002	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

* - used only with TURIOS LE + SE19

				-
Armrest SE53 Turios				Plan No.
				EP00488/J
Scale: 1:1	Drafted	Verified	Confirmed	
Date	12.02.2020	7.11.2025	7.11.2025	
Initials	DR	FM	MK	

Electrical components presented in 2992521/2992522/2992523/2992524/2992525