



#	Reference	AS	Description	Option	Qty
1	EC_BAT1, EC_BAT2	1917939	EC BAT BS E BATTERIJ 12V DAB12-44EV-HD AGM	44AH	2
		1916730	EC BAT 50AH CPX	PL	2
		1904684	EC BAT 50AH MK	50AH	2
2	EC_A2	2995296	EC A TURIOS BATTERY HARNESS SET 1		1
3	EC_A1	2995225	EC A TURIOS BATTERY FUSE 150A		1
4	EC_COM1	1910896	EC COM DLS-PM50-A	6Km/h	1
		1909020	EC COM DLX-REM050-A	REM050	1
		1910700	EC COM DLX-REM060-A	REM060	1
6	EC_MOT1	1910945	EC MOT 6Km/h R	6Km/h	1
7	EC_MOT2	1910944	EC MOT 6Km/h L	6Km/h	1
8	EC_WIRE1	1909022	EC COM GLM-BUS100-A		1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER PROCESSES WHICH WRITEN AUTHORIZATION FROM VERMEIREN GROUP

2992521/B

EC A TURIOS LE

Plan No.

2992521/B

Size: A4	Drafted	Verified	Confirmed
Date:	2025-12-01	2025-12-01	2025-12-01
Initials:	FM	MK	PM



BK BLACK
RD RED
LBL LIGHT BLUE
BL BLUE
LGR LIGHT GREEN
GR GREEN
VI VILLET
OR ORANGE
YE YELLOW
WH WHITE
LRD LIGHT RED
PK PINK
GY GRAY
/ AND
I OR