


9 EC\_A3


5 EC\_COM2

8 EC\_WIRE1  
Total length: 1000



7 EC\_MOT2

6 EC\_MOT1




#	Reference	AS	Description	Option	Qty
1	EC_BAT1, EC_BAT2	1917939	EC BAT BS E BATTERYJ 12V DAB12-44EV-HD AGM	44AH	2
		1916730	EC BAT 50AH CPX	PL	2
		1904684	EC BAT 50AH MK	50AH	2
2	EC_A2	2995296	EC A TURIOS BATTERY HARNESS SET 1		1
3	EC_A1	2995225	EC A TURIOS BATTERY FUSE 150A		1
4	EC_COM1	1918247	EC COM DLX-PM60-C	6Km/h	1
		1912302	EC COM DLX-PM75AL-C	10Km/h	1
5	EC_COM2	1912301	EC COM DLX-REM210-A	REM210	1
		1918980	EC COM DLX-REM150-A	REM150	1
		1910945	EC_MOT 6Km/h R	6Km/h	1
6	EC_MOT1	1916235	EC_MOT 10Km/h R	10Km/h	1
		1910944	EC_MOT 6Km/h L	6Km/h	1
		1916234	EC_MOT 10Km/h L	10Km/h	1
8	EC_WIRE1	1909022	EC COM GLM-BUS100-A		1
		2995292	EC A CABLE LINX SE19 SE39	SE19	1
10	EC_A4	2995293	EC A TURIOS LIGHTS E CABLE	SE19	1
11	EC_A5	2995289	EC A TURIOS LIGHTS FRONT	SE19	1
		2995290	EC A TURIOS LIGHTS REAR	SE19	1
12	EC_A6	1914000	EC COM DLX-ACU200-A	SE09	1
		1916653	EC COM DLX-USB02-A	SE06	1

- BK BLACK
- RD RED
- LBL LIGHT BLUE
- BL BLUE
- LGR LIGHT GREEN
- GR GREEN
- VI VILLET
- OR ORANGE
- YE YELLOW
- WH WHITE
- LRD LIGHT RED
- PK PINK
- GY GRAY
- / AND
- I OR

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER PROVIDES WRITTEN AUTHORIZATION FROM VERMEIREN GROUP

2992523/B

**EC A TURIOS EL**

Size: A4	Drafted	Verified	Confirmed	
Date:	2025-12-02	2025-12-02	2025-12-02	
Initials:	FM	MK	PM	

Plan No.  
**2992523/B**