

EP00156 / M

DRIVING GEAR WITH NEXUS 3 2219				Plan No.
PAGE 1/2				EP00156/ M
Scale: 1:6	Drafted	Verified	Confirmed	
Date	05/06/2014	18/03/2024	18/03/2024	
Initials	AZ	PM	AZ	

Material			
No.	Plan	Article	Quantity
1	YP03014 (Fix Wheel)	2822000	1
2	-----	-----	-
3	YP04706	2886500	1
4	6x6x50 DIN6885A	GP19000	1
5	XP02362	GP17100	1
6	XP03236	2631200	1
7	XP02365	GP17200	1
8	XP01283	2992321	2
9	M10 DIN 125-1A	1001558	5
10	M10 DIN 934	1001672	6
11	D20 DIN471	1995259	4
12	5x5x40 DIN 6885A	1940800	1
13**	---	1902733	3
14	XP01296	2992350	1
15	D160mm	2992382	1
16	Wheel 24" Left	2994426	1
17	M5x16	1006112	6
18	M8x25	1990501	2
19	---	1994568	4
20	M8 DIN 9021	1996574	2
21	M10 DIN 439B	1903216	1
22	M5x10 DIN 912/8.8	1001488	1
23	Wheel 24" Right	2994425	1
24	Chain	2993196	1,66 [m]
25	M4 DIN 985	1001678	6
26	XP03143	2996710	1
27	XP01177	2992207	1
28	X002282	1992854	6
29	M4x16	2990185	6
30*	Lagoon	1902734	3
	Vintage	1907965	3


Nexus 3 subparts - 2996710 (2992142***)			
No.	Plan	Article	Description
26a	YP04756	FP02286	nexus 3 hub with extra drilled holes
26b	---	2996706	shifting lever
26c	---	2996707	components for hub
26d**	---	2996708	clamp 5/8"
26e	---	2996709	gear 18T (2,3mm)

* - option/opcja

** - not depicted

*** - non-modified set, non-usable for clients

EP00156/M

DRIVING GEAR WITH NEXUS 3 2219				Plan No.
PAGE 2/2				EP00156/ M
Scale: 1:6	Drafted	Verified	Confirmed	
Date	05/06/2014	18/03/2024	18/03/2024	
Initials	AZ	PM	AZ	