



Material			
No.	Plan/Article	Description	Quantity
1	EP00458	LOWER FRAME	1
2	EP00455	MOTOR SYSTEM	1
3	EP00459	BATTERY BOX	1
4	EP00470	BASIC BACKREST AND SEAT	1
	EP00457	ELECTRIC BACKREST AND SEAT	1
	EP00624	ELE.ADJ. SEAT SE30 OR SE40 WITH SE36	1
	EP00625	MANUAL ADJ. SEAT SE30 AND SE41	1
5	EP00450	BI5 ARMREST	1
6	EP00461	B06 FOOTREST	1
	E00055257	BZ7/8 FOOTREST	1
* 7	EP00462	TOP CASE SE27	1
* 8	EP00477	PAVEMENT CLIMBER	1
* 9	EP00487	CRUTCH HOLDER B31	1
* 10	EP00488	ARMREST SE53	1
* 11	EP00489	WIDE ARMREST SE53	1
* 12	EP00490	MIRROR SE74	1
* 13	EP00491	CUP HOLDER	1
* 14	EP00492	STICKER CONFIGURATION	1
* 15	EP00497	EXTENDED FRAME BASIC	1
* 16	EP00498	EXTENDED FRAME ELECTRIC	1
* 17	E001090	SE08	1
	E000958	SE08 LE	1
* 18	E00061546	SE64	1
**19	E00114194	HAND BRAKE	1
* 20	E00120998	SE47	1

* - Not depicted, option,

** - Option for Germany.

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

TURIOS EV

Scale:	1:10	Drafted	Verified	Confirmed
Date	27/11/2019	27/01/2026	27/01/2026	
Initials	AS	AZ	MK	

Plan No.

EP00475/I

