



#	Reference	Art number	Description	Qty
1	EC_CON1, EC_CON51	L000007	EC CON CABLE SHOE M6 4-6MM ² HOLE 6MM	2
2	EC_CON52, EC_CON57	1991692	EC CON PLUG SB 50 GREY 2P TRACER	2
3	EC_HST1	1901731	EC HST PLASTIC CAP FOR BATTERY POSITIVE RED	1
4	EC_HST2	1994180	EC HST PLASTIC CAP FOR BATTERY NEGATIVE BLACK	1
5	EC_HST3	H005006	EC HST ISO ISOLATE COVER D10 PVC BLACK	0,35M
6	EC_WIRE9	L000034	EC WIRE 6mm ² RD	0,4M
7	EC_WIRE10	L000035	EC WIRE 6mm ² BK	0,4M
8	EC_WIRE14	1994573	EC WIRE CONNECTOR KIT POWER MODULE+LIGHTS DX2	1

- BK BLACK
- RD RED
- LBL LIGHT BLUE
- BL BLUE
- LGR LIGHT GREEN
- GR GREEN
- VI VILLET
- OR ORANGE
- YE YELLOW
- WH WHITE
- LRD LIGHT RED
- PK PINK
- GY GRAY
- / AND
- | OR

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP

PP00038_A_Turios_Battery_Harness

Size: A4	Drafted	Verified	Comformed
Date:	2025-07-30	2025-09-30	2025-09-30
Initials:	FM	MK	PM

E00008220/B

Plan No.

E00008220/B