



\* - used only with TURIOS LE + SE19

EP00489 / J

SE53 WIDE ARMREST / TURIOS				Plan No.
				EP00489 / J
Scale: 1:5	Drafted	Verified	Confirmed	
Date	12/02/2020	5/11/2025	5/11/2025	
Initials	DR	FM	MK	

Material				Material			
No.	Plan	Article	Quantity	No.	Plan	Article	Quantity
1	-	2606454	1	28*	M4X12	2994126	2
2	Y003626	1476554	1	29	-	-	-
3	Y009175	1641654	3	-	-	-	-
4	Y009176	1641754	1	-	-	-	-
5	YP08065	2605954	2	-	-	-	-
6A	YP08058	2605454 (L)	1	-	-	-	-
	YP08059	2605554 (R)	1	-	-	-	-
6B*	YP08060	2605654 (L)	1	30	M5x16	1907493	2
	YP08061	2605754 (R)	1	31	-	-	-
7	-	-	-	32	Spring plunger	1991432	1
8	YP09265	2642654	1	33	-	-	-
9	-	-	-	34a	Armrest wide Soft Gray	1999993 (R)	1
10	X005565	1912287	1		Armrest wide Soft Gray	1999992 (L)	1
11	Armrest right	1913465	1		Armrest wide Soft Black	2998254 (R)	1
	Armrest left	1913464	1		Armrest wide Soft Black	2998253 (L)	1
12	XP00227	2990147	2/4**	34b	L30	1998790 (R)	1
13	-	-	-		L30	1998789 (L)	1
14	Washer M6	1000973	8	35*	Washer M4 DIN127B	1998515	2
15	M5x12	1001489	2	36	M6x45	2994058	2
16	Tube Cap 20x20	1001612	2	37	Tube Cap 16x16	2994195	2
17	Nut M10	1001672	1	38	XP03080	2996656	1
18	Nut M5	1001679	6				
19	Nut M6	1001680	6				
20	M6x40	1003111	2				
21	Star Knob	1003658	1				
22	M5x20	1003907	4				
23	M6x40	1006111	4				
24	-	-	-				
25	M6x16	1008459	2				
26	Nut M6	1008492	2				
27	-	-	-				

Bill of Material 1913465/1913464 for EP00489			
No.	Plan	Article	Quantity
11a	X005567	1912286 (R)	1
	X005566	1912285 (L)	1
11b	X005517	2014999	1
11c	X005561	2021699	1
11d	M6x16	1008459	1
11e	Nut M6	1008492	1
11f	Rivet Ø4	1992805	2
11g	M5x16	3991136	2
11h	Washer M5	1906561	2
11i	Sticker	T001002	1

\* - used only with TURIOS LE + SE19

\*\* - for L30

Electrical components presented in 2992521/2992522/2992523/2992524/2992525

EP00489 / J

SE53 WIDE ARMREST / TURIOS				Plan No.
				EP00489 / J
Scale: 1:5	Drafted	Verified	Confirmed	
Date	12/02/2020	5/11/2025	5/11/2025	
Initials	DR	FM	MK	