



#	Reference	Art number	Description	Qty
1	EC_WIRE10, EC_WIRE14	1907749	EC WIRE CABLE 4X0,5 BLACK	1,6M
2	EC_L1, EC_L3	1907340	EC L LED HEAD LIGHT	2
3	EC_L2, EC_L4	1991679	EC L LIGHT DIRECTION INDICATOR	2
4	EC_CON47-EC_CON50	2990811	CABLE CHOE 68/76 MAIL	4
5	EC_WIRE8, EC_WIRE9	L000022	EC WIRE 0.75MM ² CABLE GR/YE	0,1M
6	EC_CON39, EC_CON41, EC_CON42, EC_CON45	2996202	CABLE SHOE 68/76 FEMAL	4
7	EC_CON40, EC_CON43, EC_CON44, EC_CON46	L000079	EC CON CABLE SHOE 6.3 1.5-2.5MM ² BLUE ISOLATED	4
8	UNDEPICTED	L000081	CABLE SHOE 6.3 1.5-2.5MM ² BLUE ISOLATED	2
9	UNDEPICTED	W620150	HANDY BUNDLER 0.04M BLACK	0,1M
10	UNDEPICTED	W620151	CABLE BINDER SEAL BLACK	2

BK BLACK
 RD RED
 LBL LIGHT BLUE
 BL BLUE
 LGR LIGHT GREEN
 GR GREEN
 VI VILLET
 OR ORANGE
 YE YELLOW
 WH WHITE
 LRD LIGHT RED
 PK PINK
 GY GRAY
 / AND
 I OR

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS
 PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN
 WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE,
 EXCEPT WHERE SHOWN USER PROCESSOR'S WRITTEN AUTHORIZATION
 FROM VERMEIREN GROUP

2995290/B

Plan No.

2995290/B

EC A LIGHT LED REAR LINX TURIOS

Size: A4	Drafted	Verified	Conformed
Date:	2025-07-30	2025-09-30	2025-09-30
Initials:	FM	MK	PM

