



#	Reference	Art number	Description	Qty
1	EC_WIRE11, EC_WIRE12	L000068	EC WIRE 6mm2 BR	0,5M
2	EC_FUS2	1910488	EC FUS Fuse 150A	1
3	EC_HST15, EC_HST17	1994180	EC HST PLASTIC CAP FOR BATTERY NEGATIVE BLACK	2
4	EC_HST16, EC_HST21	1901731	EC HST PLASTIC CAP FOR BATTERY POSITIVE RED	2
5	UNDEPICTED	1007000	HEX SOCKET BUTTON HEAD SCREW M6X12	2
6	UNDEPICTED	1004414	LOCKWASHER DIN127 M6	2
7	UNDEPICTED	1001680	SELFLOCKING NUT DIN985 M6 ZINC	2
8	EC_CON21, EC_CON22	L000007	EC CON CABLE SHOE M6 4-6MM ² HOLE 6MM	2
9	EC_CON18-EC_CON20, EC_CON23	L000001	EC CON CABLE SHOE M6 4-6MM ²	4
10	EC_HST18, EC_HST20	T001009	EC HST Heat Shrink Tubing D19	0,1M

BK	BLACK
RD	RED
LBL	LIGHT BLUE
BL	BLUE
LGR	LIGHT GREEN
GR	GREEN
VI	VILLET
OR	ORANGE
YE	YELLOW
WH	WHITE
LRD	LIGHT RED
PK	PINK
GY	GRAY
/	AND
	OR

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORIZATION FROM VERMEIREN GROUP.

EC A TURIOS BATTERY FUSE 150A

Size: A4

Drafted

Verified

Comformed

Date:

2025-07-30

2025-09-30

2025-09-30

Initials:

FM

MK

PM

2995225/A

Plan No.

2995225/A

