

Sheet 1/2

2991818

E000857 / E


| | | | | |
|---------------------------------------|------------|------------|------------|-------------|
| <h3>TIMIX SEAT FRAME GT ELECTRIC</h3> | | | | Plan No. |
| | | | | E000857 / E |
| Scale: 3:20 | Drafted | Verified | Confirmed | |
| Date | 25/01/2013 | 19/09/2016 | 19/09/2016 | |
| Initials | SD | SD | PJ | |

| Material | | | | | | | |
|----------|-----------|--------------|----------|-----|----------------|--------------|----------|
| No. | Plan | Article | Quantity | No. | Plan | Article | Quantity |
| 1 | X000862 | 1006463 | 4 | 14 | M8 DIN125 | 1001556 | 6 |
| 2 | X001016 | 1001682 | 10 | 15 | M6 DIN985 | 1001680 | 12 |
| 3 | M12 PA6.6 | 1993362 | 2 | 16 | M10 DIN985 | 1001683 | 2 |
| 4 | X003403 | 1617400 | 2 | 17 | M6x12 ISO7380 | 1007000 | 4 |
| 5 | X004860 | 1906817 (40) | 1 | 18 | M10x35 ISO7380 | 1009442 | 1 |
| | X004855 | 1906741 (45) | 1 | 19 | M8x50 ISO7380 | 1902483 | 2 |
| | X004887 | 1906818 (50) | 1 | 20 | M8x20 DIN7991 | 1903136 | 4 |
| 6 | X004888 | 1906839 (40) | 1 | 21 | --- | 1906790 | 1 |
| | X004857 | 1906740 (45) | 1 | 22 | --- | 1906843 (40) | 1 |
| | X004889 | 1906841 (50) | 1 | | --- | 1906844 (45) | 1 |
| 7 | Y017988 | 1774154 | 1 | | --- | 1906845 (50) | 1 |
| 8 | Y018940 | 1778754 (40) | 1 | 23 | --- | 1990648 | 2 |
| | Y018164 | 1776654 (45) | 1 | 24 | --- | 1991481 | 1 |
| | Y018361 | 1778654 (50) | 1 | 25 | M8x16 DIN912 | 1001498 | 4 |
| 9 | Y018564 | 1774200 | 2 | 26 | --- | 2993352 | 1 |
| 10 | Y017414 | 1775254 (40) | 1 | 27 | M10x45 ISO7380 | 1908848 | 1 |
| | Y017984 | 1769754 (45) | 1 | 28 | M8x30 ISO7380 | 1992717 | 2 |
| | Y017889 | 1776054 (50) | 1 | 29 | M6x30 DIN7984 | 1992721 | 8 |
| 11 | Y018411 | 1778354 (40) | 1 | 30 | --- | 1997922 | 2 |
| | Y018368 | 1778454 (45) | 1 | 31 | --- | 2990413 | 2x0,4 m |
| | Y018309 | 1778554 (50) | 1 | 32 | --- | 1906791 | 1 |
| 12 | Y019511 | 1784654 (40) | 1 | 33 | --- | 1906792 | 1 |
| | Y019510 | 1784754 (45) | 1 | 34 | M8x20 ISO7380 | 1992349 | 2 |
| | Y019512 | 1784854 (50) | 1 | 35 | M8 DIN125 | 1003199 | 2 |
| 13 | Y019513 | 1784954 (40) | 1 | | | | |
| | Y019509 | 1785054 (45) | 1 | | | | |
| | Y019514 | 1785154 (50) | 1 | | | | |

Sheet 2/2

2991818

E000857 / E

| | | | | |
|------------------------------|------------|------------|------------|---|
| TIMIX SEAT FRAME GT ELECTRIC | | | | Plan No. |
| | | | | E000857 / E |
| Scale: 3:20 | Drafted | Verified | Confirmed |  |
| Date | 25/01/2013 | 19/09/2016 | 19/09/2016 | |
| Initials | SD | SD | PJ | |