



BK	BLACK
RD	RED
LBL	LIGHT BLUE
BL	BLUE
LGR	LIGHT GREEN
GR	GREEN
VI	VILLET
OR	ORANGE
YE	YELLOW
WH	WHITE
LRD	LIGHT RED
PK	PINK
GY	GRAY
/	AND
	OR

#	Reference	AS	Description	Option	Qty
1	EC_BAT1, EC_BAT2	1917939	EC BAT BS E BATTERIJ 12V DAB12-44EV-HD AGM	44AH	2
		1916730	EC BAT 50AH CPX	PL	2
		1904684	EC BAT 50AH MK	50AH	2
2	EC_A2	2995297	EC A TURIOS BATTERY HARNESS SE19	SE19	1
3	EC_A1	2995225	EC A TURIOS BATTERY FUSE 150A		1
4	EC_COM1	1910896	EC COM DLS-PM50-A	6Km/h	1
5	EC_COM2	1909020	EC COM DLX-REM050-A	REM050	1
1910700		EC COM DLX-REM060-A	REM060	1	
6	EC_MOT1	1910945	EC_MOT 6Km/h R	6Km/h	1
7	EC_MOT2	1910944	EC_MOT 6Km/h L	6Km/h	1
8	EC_WIRE1	1909022	EC COM GLM-BUS100-A		1
9	EC_A3	2995295	EC A TURIOS LIGHT CONTROLLER	SE19	1
10	EC_A4	2995293	EC A TURIOS LIGHTS E CABLE	SE19	1
11	EC_A5	2995289	EC A TURIOS LIGHTS FRONT	SE19	1
12	EC_A6	2995290	EC A TURIOS LIGHTS REAR	SE19	1

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER PROVIDES WRITTEN AUTHORIZATION FROM VERMEIREN GROUP.

2992522/B

EC A TURIOS LE + SE19

Plan No.

2992522/B

Size: A4	Drafted	Verified	Confirmed
Date:	2025-12-02	2025-12-02	2025-12-02
Initials:	FM	MK	PM

