



Material			
No.	Plan	Article	Quantity
1	E001038	1913465 (R)	1
		1913464 (L)	1
2	P00046914	2997543	1
3	M6x55 DIN912	1908726	2
4	M6	1008492	2
5	M5X20 DIN912	1003907	4
6	M5 DIN985	1001679	4
7	XP03080	2996656	1
8	M6 DIN125A	1001564	2

Method PL		
Procedure	Description	Tool
Drill	2 holes Ø6	500-052

Only left side depicted
Quantities for one side

BI5 ARMREST ECLIPS X4				2996746
				Plan No. E00046972/A1.1
Scale: 3:10	Drafted	Verified	Confirmed	
Date	2022-11-04	04/11/2022	04/11/2022	
Initials	AS	AZ	AS	