



#	Reference	Art number	Description	Qty
1	EC_WIRE7	1009355	EC WIRE CABLE 4X0,75 BLACK	0,35M
2	EC_CON15	1912300	EC CON PLUG SET LINX SE19	1
3	EC_CON11-EC_CON14	2996202	CABLE SHOE 68/76 FEMALE	4
4	EC_HST11	T001009	EC HST Heat Shrink Tubing 19	0,4M

BK	BLACK
RD	RED
LBL	LIGHT BLUE
BL	BLUE
LGR	LIGHT GREEN
GR	GREEN
VI	VILLET
OR	ORANGE
YE	YELLOW
WH	WHITE
LRD	LIGHT RED
PK	PINK
GY	GRAY
/	AND
	OR

NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORIZATION FROM VERMEIREN GROUP.

2995292/B

EC A LINX SE19 SE39 E CABLE

Plan No.

2995292/B

Size: A4	Drafted	Verified	Comformed
Date:	2025-07-30	2025-09-30	2025-09-30
Initials:	FM	MK	PM

