



A	SEAT DEPTH	380, 400, 420, 440, 460, 480, 500
B	SEAT WIDTH	390, 420, 450, 480, 510
C	SEAT ANGLE	0°, 2°, 4°, 6°
D	FOOTREST LENGHT	347, 389, 424, 449, 474,
E	ARMREST HEIGHT	245, 260, 275, 290, 305, 320
F	WIDTH OF ARMREST	411, 471, 531 width adjustment
G	SEAT HEIGHT	440, 490, 540
H	FOOTPLATE HEIGHT	180
I	BACKREST ADJ	0°, 8°, 16°, 24°, 32°



NOTICE: THE INFORMATION CONTAINED ON THIS DRAWING IS PROPRIETARY TO VERMEIREN GROUP AND SHALL NOT BE DISCLOSED IN WHOLE OR IN PART, OR USED FOR ANY DESIGN OR MANUFACTURE, EXCEPT WHEN SUCH USER POSSESSES DIRECT WRITTEN AUTHORISATION FROM VERMEIREN GROUP.

DIMENSIONS TURIOS BASIC				---
				Plan No. EP00483/C
Scale: 1:10	Drafted	Verified	Confirmed	
Date	16/01/2020	18/08/2025	18/08/2025	
Initials	AS	AZ	MK	