

Eclips X4 30°

30° reclinable backrest

Lightweight aluminum wheelchair with a reclining backrest of 30°, reinforced cross and Quickrelease wheels. Standard equipped with short armrests for easy table access.



VERMEIREN
we care for you

Standard

- Height adjustable handles
- Inclination of the back by hydraulic springs for more comfort during the tilting
- Foldable and removable armrests
- Rear wheels with Quickrelease system
- Internal and external collapsible footplates

Colours



Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.

Quality guaranteed

Regarding quality, the **Eclips X4 30°** complies with the highest norms, laid down by certified test institutes in most European countries.



Different adjustment possibilities for the seat height and seat depth



Adjustment of the backrest by means of gas spring for a smooth adjustment

Technical data

■ Total width	590/620/640/660/680/700/740 mm
■ Total height	910 - 1015 mm
■ Total length	1000 - 1210 mm
■ Folded width	300 mm
■ Seat width	390/420/440/460/480/500/550 mm
■ Height armrests	220 mm

■ Seat height	390 - 520 mm
■ Height backrest	420 - 460 mm
■ Seat depth	420 - 500 mm
■ Weight	19,5 kg
■ Maximum user weight	150 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you