



Material			
No.	Plan	Article	Quantity
1	Y003413	1463254	1
2	Y003414	1463354	1
3	Y003416	B000272	1
4	Y007779	1607904	1
5	VGf M8x15	1006661	2
6	15x15	1005899	1
7	M6x20	1006999	4
8	M6x12	1007000	4
9	VGf40 M8x16	2990024	1
10	Y017183	1748154	1
11*	Soft Black	1919072	1
	Soft Grey	1907019	1
	Cordura	1912437	1
	Black Skai	1911716	1

*not depicted

1907195

E000987 / F

L55 TYPE INOVYS II / V-SERIES				Plan No.
				E000987 / F
Scale: 1:3	Drafted	Verified	Confirmed	
Date	18/06/2015	14/12/2023	14/12/2023	
Initials	TT	PM	AS	